

T3 Uptake	33	24-39 %	01
Free Thyroxine Index	1.5	1.2-4.9	01
Testosterone, Free/Tot Equilib			
Testosterone, Serum	209	LOW 264-916 ng/dL	01
Adult male reference interval is based on a population of healthy nonobese males (BMI <30) between 19 and 39 years old.			
Travison, et.al. JCEM 2017,102;1161-1173. PMID: 28324103.			
Testosterone, Free	8.07	5.00-21.00 ng/dL	01
% Free Testosterone	3.86	1.50-4.20 %	01
Prostate-Specific Ag, Serum			
Prostate Specific Ag, Serum	0.5	0.0-4.0 ng/mL	01
Roche ECLIA methodology.			
According to the American Urological Association, Serum PSA should decrease and remain at undetectable levels after radical prostatectomy. The AUA defines biochemical recurrence as an initial PSA value 0.2 ng/mL or greater followed by a subsequent confirmatory PSA value 0.2 ng/mL or greater.			
Values obtained with different assay methods or kits cannot be used interchangeably. Results cannot be interpreted as absolute evidence of the presence or absence of malignant disease.			
IGF-1			
Insulin-Like Growth Factor I	190	75-216 ng/mL	01
Estradiol, Sensitive			
Estradiol, Sensitive	18.8	8.0-35.0 pg/mL	01
This test was developed and its performance characteristics determined by LabCorp. It has not been cleared by the Food and Drug Administration.			
Methodology: Liquid chromatography tandem mass spectrometry (LC/MS/MS)			

Performing Laboratory Information:

01: LabCorp Burlington, 1447 York Court, Burlington NC, phone: 800-762-4344
 Medical Director: MD William F Hancock